



Phoniro Release Information

When we update our IT system Phoniro Care and our apps, you can find what we have done in the Release Information. We continuously develop our applications to give you better functions and higher stability. By keeping up-to-date, you will get the most from your investment.

When we launch a new version of Phoniro Care, any of our apps, and/or software for our alarm servers, it is called a release. These new versions we collect into a release package. Sometimes there are several new versions (releases) in one release package but other times there is only one new version included.

Versions included in this release package

Phoniro Care version

Phoniro Care
3.61 (current)

Alarm server versions

Phoniro Local Alarm Server (PLAS)	Phoniro Alarm Server Cloud
1.15.2 (current)	1.15.2 (new release)

Apps versions

Phoniro Home Care (Android & iOS)	Phoniro PI (Android)	Senior Living (Android)
3.0 (new release)	3.0 (new release)	2.1 (current)

What features are included in a release?

Some of the features will be automatically shown to you after an upgrade.

However, other features need to be activated and configured for them to take effect. Such features are marked with this symbol: ⚙️. These features might involve configuration costs as well as a monthly cost.

Contact our Service Center for more information on how the features can be used in your organization.


Digital Key Management

Updated platform for our apps, Home Care and Phoniro PI


We have updated the platform on which the Home Care app is built. This should not result in any changes to the app's functionality, but it may cause minor visual changes compared to the previous version. For example, the text size and the size or shape of buttons may have changed

Ability to choose text size


The text size is slightly smaller compared to previous versions of Home Care and PI. However, users now have the option to adjust the text size in the app themselves under **Main Menu > Settings > Font size**.



View in Home Care to reach Font size settings



View in Phoniro PI to reach Font size settings



Font size settings in both apps

How do I get access to new releases?

Phoniro Care

For customers with Phoniro Care in their own operation

- You can download the files you need, together with technical information, from the storage space you have gained access to from us. Then you upgrade on your own.

For customers who buy Phoniro Care as a service

- We upgrade Phoniro Care for you and make sure you have the latest software. Easy!

Do you know that we offer fast and free upgrades when you let us take care of the operation of Phoniro Care? Contact us for more information.

The apps Home Care, Phoniro PI and Senior Living

- Our apps are usually continuously updated via Google Play/AppStore.

Would you like to learn more?

More detailed information about the features is provided in our user manuals. All Phoniro Care's solutions are described on our website <https://www.phoniro.com/>